



# STERILON










## Monofilament Polyamide

### Non Absorbable Suture USP












<b>Made of</b>	Polyamide 6.6
<b>Color</b>	Black or Blue
<b>Absorption</b>	None
<b>Characteristics</b>	Characterized by high strength and good knotting abilities, smooth surface does not support bacterial growth
<b>Tensile strength</b>	High
<b>Knotting strength</b>	Good
<b>Tissue reaction</b>	Minimal
<b>Sterilization method</b>	Ethylene oxide or Irradiation
<b>Packaging</b>	One, Two & Three dozen wise, in environmental friendly cardboard boxes
<b>Shelf life</b>	5 years
<b>CE-mark</b>	CE 0434
<b>Standard sizes</b>	USP 10/0 to 2

Sterilon sutures listed in this catalogue conform to the specifications of the United States Pharmacopeia .






## STERILON (Needled Sutures)

Needle Length	Description	Suture Length	Size			
			10-0 (0.2)	9-0 (0.3)	8-0 (0.4)	
<b>6mm</b>						
 	3/8 Circle Spatulated (Premium Point)	38 cm	SFN 3718			
 	3/8 Circle Spatulated (Premium Point)	38 cm		SFN 3715	SFN 3322	
DOUBLE ARMED   	3/8 Circle Spatulated (Premium Point)	38 cm	SFN 3719	SFN 3714		
DOUBLE ARMED   	3/8 Circle Advanced Spatulated (Premium Point)	38 cm	SFN 3717			
<b>8mm</b>						
 	3/8 Circle Reverse Cutting (Premium Point)	38 cm			SFN 3343	

# STERILON (Needled Sutures)

Needle Length	Description	Suture Length	Size				
			6-0 (0-7)	5-0 (1)	4-0 (1.5)	3-0 (2)	
<b>10mm</b>							
▼	DS 	3/8 Circle Reverse Cutting	38 cm			SFN 3326	
			70 cm	SFN 3320	SFN 3323		
<b>12mm</b>							
▽	DS 	3/8 Circle Reverse Cutting (Premium Point)	70 cm		SFN 3317		
▲	DC 	3/8 Circle Cutting	35 cm		SFN 3315		
▽	DS 	3/8 Circle Reverse Cutting (Premium Point)	70 cm			SFN 3318 SFN 3321	
<b>16mm</b>							
●	DR 	3/8 Circle Round Bodied	70 cm		SFN 3424		
<b>17mm</b>							
▼	DS 	3/8 Circle Reverse Cutting	70 cm			SFN 3931	
<b>26mm</b>							
▲	DC 	3/8 Circle Cutting	70 cm	SFN 3331			
▼	DS 	3/8 Circle Reverse Cutting	70 cm			SFN 3328	
Needle Length	Description	Suture Length	Size				
			3-0 (2)	2-0 (3)	1-0 or 0 (3.5)	1 (4)	2 (5)
<b>40mm</b>							
LOOP	HR 	1/2 Circle Round Bodied	150 cm			SFN 3340	
●							
HEAVY	HR 	1/2 Circle Round Bodied	100 cm			SFN 3346	
●							SFN 3347
<b>45mm</b>							
▼	DS 	3/8 Circle Reverse Cutting	70 cm		SFN 3336		

# STERILON (Needled Sutures)

Needle Length	Description	Suture Length	Size				
			3-0 (2)	2-0 (3)	1-0 or 0 (3.5)	1 (4)	2 (5)
<b>45mm</b> HEAVY 	1/2 Circle Reverse Cutting	100 cm				SFN 3338	
<b>50mm</b> HEAVY LOOP 	1/2 Circle Round Bodied	150 cm				SFN 3348	
<b>60mm</b> 	Straight Cutting	70 cm	SFN 3388	SFN 3390	SFN 3389		
<b>90mm</b> 	3/8 Circle Cutting	70 cm				SFN 3397	SFN 3398
	3/8 Circle Cutting	75 cm					SFN 3398A

# STERILON (Non-Needled Sutures)

Length	Packaging	Size (U.S.P.)				
		3-0 (2)	2-0 (3)	1-0 or 0 (3.5)	1 (4)	2 (5)
76 X 2 cm	Sterile Foil	SMP 712	SMP 713	SMP 714	SMP 71	SMP 716

\* Needle to Suture Combinations are available as per customer specifications.