

# STERICRYL


*Monofilament Synthetic Absorbable Suture EP*  
(Poliglecaprone 25)

# STERICRYL











## Monofilament Poliglecaprone 25 Synthetic Absorbable Suture

<b>Made of</b>	75% Polyglycolide-co-Polycaprolactone 25%
<b>Color</b>	Violet or Undyed
<b>Material Type</b>	Synthetic , Monofilament
<b>Absorption</b>	Via Hydrolysis
<b>Characteristics</b>	<b>Better substitute to Catgut.</b> Unprecedented pliability and smooth passage through tissues, easy handling & knotting. <b>Dyed</b> suture provides wound support upto 28 days and <b>Undyed</b> upto 21 days
<b>Tensile strength</b>	VIOLET : 60-65% at 7 days, 25-30% at 14 days, UNDYED : 50-55% at 7 days, 20-25% at 14 days
<b>Knotting strength</b>	Good
<b>Tissue reaction</b>	Minimal
<b>Sterilization method</b>	Ethylene oxide
<b>Packaging</b>	One, Two & Three dozen wise, in environmental friendly cardboard boxes
<b>Shelf life</b>	5 years
<b>CE-mark</b>	CE 0434
<b>Standard sizes</b>	EP 6/0 to 1

Stericryl Sutures listed in this catalogue conform to the specifications of the European Pharmacopeia.



# STERICRYL (Needled Sutures)

Needle Length	Description	Suture Length	Size			
			4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 or 0 (3.5)
<b>16mm</b>						
	3/8 Circle Reverse Cutting (Premium Point)	70 cm Undyed	STR 1205			
	3/8 Circle Round Body	70 cm Dyed Violet	STR 1748			
<b>20mm</b>						
	1/2 Circle Round Body	70 cm Dyed Violet		STR 1737		
<b>25mm</b>						
	3/8 Circle Cutting	70 cm Dyed Violet		STR 1726		
	3/8 Circle Cutting (Premium Point)	70 cm Undyed		STR 1326		
<b>26mm</b>						
	1/2 Circle Oval Round Body	70 cm Undyed		STR 1648		
	1/2 Circle Oval Round Body (JB Needle)	70 cm Undyed		STR 1664		
	1/2 Circle Oval Round Body	70 cm Undyed			STR 3665	
	1/2 Circle Oval Round Body (JB Needle)	70 cm Dyed Violet		STR 1764		
<b>30mm</b>						
	1/2 Circle Round Body	70 cm Dyed Violet			STR 1741	STR 1742

\* Needle to Suture Combinations are available as per customer specifications.